



# *Advanced Personal Laser Hair Remover*

USER MANUAL

DLX-1088



This product **MUST** be used responsibly.

Please read **ALL** instructions including the User Manual **CAREFULLY** before use.



CLASS 1  
LASER PRODUCT

## Package Includes

1. Laser Hair Removal System Unit x 1
2. Switching Power Adapter x 1
3. English Instruction Manual x 1

## Warning

Please read all warnings and safety precautions before use, and strictly follow them when using DLX-1088 Advance Personal Laser Hair Remover. DLX-1088 is a personal use device intended for the removal of unwanted hair. DLX-1088 can be used by individuals with light skin tones and brown-black hairs.

Before using DLX-1088 for the first time, please read this User Manual in its entirety. Pay special attention to each section on device use procedures, device operation, and after-user procedures. We recommend you re-check and re-understand yourself with this User Manual before each use of DLX-1088.

# Skin Tone Chart

Check your skin tone : Using the chart on the chart below. Find the corrected row of the chart that best matches with your skin in the area you want to treat. Don't use the DLX-1088 if your skin is too dark.

## Intended for These Skin Tones

White
Ivory
Beige
Light Brown

## Not Safe for These Skin Tones - DO NOT USE

Medium Brown
Dark Brown/Black

## Table of contents

1. What is DLX-1088 Advance Personal Laser Hair Remover .....	4
2. Advantages .....	4
3. Warnings & Safety information .....	4
4. Precautions .....	5
5. Warnings .....	6
6. Risks .....	6
7. Package Contents and Device Description .....	7
8. Quick guide .....	8
9. What to expect – a timeline .....	9
10. Plan your DLX-1088 for best results .....	10
11. Q & A .....	11
12. Troubleshooting .....	11
13. Taking care of your DLX-1088 system .....	12
14. Performance environment .....	12
15. Specifications .....	12
16. Customer Service .....	12

Thank you choosing DLX-1088. Now you can permanently reduce unwanted hair in the privacy of your own home using the same procedure as a beauty shop, but at a fraction of the price. The Diode Laser is the world's first affordable home laser unit and is available right now.

## 1. What is DLX-1088 Advance Personal Laser Hair Remover

The DLX-1088 system disables the follicles in your skin that produce hair. It emits a pulse of laser light, which burns up the dark pigment inside the hair. This heat disables the hair follicle. In other words the laser zaps only what's dark-that's why it works only on naturally brown or black hair.

To be effective, the laser needs to reach each hair follicle during its growth cycle. This requires treatment twice a month for 3 months and then once a month for 3-5 months, which results in gradual reduction after the third treatment.

## 2. Advantages

When used directly, the DLX-1088 system gives you long term removal of unwanted hair without the need for daily shaving, waxing, or plucking.


The results may vary from user to user and it depends on any body area you treat and how well you follow the instructions in this user manual, but the clinical study demonstrated the following advantages.

It is a much faster and easier way for a hair removal than other existing methods. Almost no painful and shorter time-consuming compared with other treatments. DLX-1088 removes your hair fast and clean.

You can feel the soft skin resulted from pore constriction after treated. Unlike shaving, it does not harm your skin and cause any inflammation. This can help minimize pores.

If you have tanned too much or have dyed the hair, wait until the skin is returned into normal color because laser has a characteristic feature that reacts with dark hair easily.

## 3. Warnings & Safety information



1. As with any electrical device, certain precautions must be taken in order to ensure your safety when using DLX-1088. DLX-1088 is an electrical device. As such it should always be kept away from water.
2. Do not use this appliance near bath-tubs, showers, basins or other vessels containing water. Don't touch the device with wet hands.

3. Keep this device out of the reach of children. The DLX-1088 should not be used on or by children below the age of 15.
4. Don't look with or without protection (optical lenses) directly to the laser beam aperture.
5. Don't attempt to open or repair your DLX-1088 device. An authorized repair center technicians must perform all repairs. Contact customer service to arrange for repair of a broken or damaged device.
6. Send the device to your dealer or service center for maintenance each year.
7. If any failure problems, please contact your dealer or service center.
8. Don't use the device during a thunderstorm and keep flammable materials away from the device.
9. Protect the device from sunlight directly.
10. Be aware of fire, an electric shock or injury.  
220 – 240Vac 50Hz 0.3A

#### 4. Precautions

1. While making the treatment, a flash red light may be seen coming from the laser applicator. Do not put your face near or expose your eyes to the flash.
2. Don't use on parts of your skin that are darker than your usual skin tone, tattoos, dark brown or black spots (such as large freckles, birthmarks, or moles), or skin that is red or swollen from a previous treatment. Any dark skin may absorb too much laser light, which can injure the skin.
3. Use a cool towel to keep your overall skin temperature down so that you can ensure the heat created by the laser for a longer period of time. It is recommended to use a Cooling gel after treatment.
4. Don't use the DLX-1088 on or around the eyes, eyebrows, or eye lashes, to prevent causes or any serious eye injury.
5. Your skin may become red right after treatment or within 24 hours. This redness generally clears up within a few hours.

## 5. Warnings

1. Don't use anywhere on the nipples, genitals or around the anus. These areas may have more sensitive skin, a darker skin tone, and/or a greater hair density, and using the DLX-1088 system may injure the skin.

### Don't use the DLX-1088 system.

2. If you have a history of skin cancer or areas of potential skin malignancies.
3. If you have allergies.
4. If you are pregnant or breast-feeding. The DLX-1088 system has not been tested on pregnant women.
5. After drinking, taking medicine, suffering from fatigue, or suffer from a high fever.
6. If you have too much sensitive skin conditions. Don't use the DLX-1088 system for any purpose other than hair removal.
7. Don't use if you are under the age of 15 or use on anyone else under the age of 18. The use of the DLX-1088 system on or by children has not been studied and could result in serious injury.
8. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
9. Children should be supervised to ensure that they do not play with the appliance.
10. The product is not working properly or problems caused by electricity static, the product will be reset automatically, and this is the normal protective mechanism in our original product design. Only press the RESET and the product can be worked normally.

## 6. Risks

1. Your skin may become red right after treatment or within 24 hours. This redness generally clears up within a few hours. See your doctor if the redness doesn't go away within 2-3 days.
2. In rare cases, some users with a particular genetic make-up may see an increase in the number of hairs in the treated area after the initial treatment.
3. If this increase persists over several treatments, stop using the DLX-1088 and see your doctor.
4. Don't use the device if the skin color changes or on areas such as blisters.


## 7. Package Content and Device Description

While open the DLX-1088 package, you will find the following parts, and please read & understand the part names before operating this product.



## 8. Quick guide


1. Clean the skin area to be treated, remove all traces of cosmetics and creams.
2. For best results, shave or clip hair to be removed to 1-2mm length. It can help laser burns the root of the hair follicle correctly.
3. Make sure the treated area has to be wiped with wet towel swipe or apply cooling gel 2 -3 minutes prior the treatment.
4. Plug in the power adapter or put 4 x AA alkaline batteries (NOT INCLUDED) and switch on the DLX-1088 home Personal Laser Hair Remover.



5. Unlock the laser with the safety key as following position.



6. The laser will be default to Lo but you can change to Med or Hi power setting by pressing the Mode Button.
7. Set the laser on Lo setting for your first treatment.
8. Hold the Handle and move the Skin Contact Bars to the skin hairs in the selected area.



9. A gently beep sound comes out and Contact LED On.
10. Press the Laser Release button, laser generates to burn hairs.
11. Pick up and move the Skin Contact Bars to another selected area for next treatment.
12. The DLX-1088 laser hair remover works best if you treat one hair at a time.
13. It is normal to feel a slight tingling sensation.



14. When treating on facial area, avoid pointing laser at the eyes.
15. After treatment, apply Cooling gel to treated areas. It will provide an instant cooling effect and proper hydration of the epidermis.
16. Don't treat the same area of skin more than once every 2 weeks. It can take up to 2 months for a hair follicle to die and release its hair.
17. Let the hairs gradually fall out, may take up to 10 to 14 days.
18. It can take up to 5 treatments to remove hairs permanently.
19. Must be followed above procedure to use the device, otherwise it will be created a Dangerous Harmful radiation.
20. Do not put your face near or expose your eyes to the Laser Aperture while using the device.
21. If any interference causes no function, press the RESET button.

## 9. What to expect – a timeline

### Before treatment

Fully wash skin area to be treated, removing all traces of cosmetics and creams.

### During treatment

Don't wax or grip hairs out. For best results, shave or clip hair to be removed to 1-2 mm length. Some users experience a warm sensation, while others experience no sensation at all. As the amount of hair in the treated areas decreases, the sensation will also decrease. Treatment sessions with DLX-1088 should be spaced every two weeks

### After the first treatment

There is no any visible changes to the hair or the appearance of the skin. Usually, some users experience mild redness that quickly dissipates within 24 hours.

### After the second month of treatment (2 – 3 treatments)

Hair growth will begin to reduce and become less dense. In some cases, it will appear as if the hair is continuing to grow. If you notice this, don't panic – this is simply the skin working to eject the treated hair from the skin.

### After three months (3 – 4 treatments)


You will begin to see a visible reduction in the amount of hairs re-growth, especially in areas where hair growth tends to be dense, like the bikini line for women or the back or chest for men.

### After six to eight months

Enjoy skin that looks and feels as though the hair were never there. Do regularly touch-ups only if necessary.

## 10. Plan your DLX-1088 for best results

Don't wax or grip hairs out. For best results, shave or clip hair to be removed to 1 - 2mm length




Do not depilate your hair, or the laser will not be effective.

Use the laser on the LOW setting for the first treatment to ascertain how your skin will react.

To get the best results from laser hair removal, use on the High setting.

Point the laser directly at the hair for permanent hair removal. Although the laser is transmitted into the hair root through the skin, it penetrates quite deeply in order to destroy it.



You may experience a slight burning or tingling sensation in use - this is normal and means that the hair root has been successfully destroyed.

Adjust the intensity of the laser depending on your sensitivity to the stinging sensation.

After 30 minutes of use, we recommend a 5 minute break to protect both your skin and the laser.

The key to long-lasting results is to treat the same area every 2 weeks over a period of 6-8 months.

For difficult-to-treat areas such as armpits, we recommend that you get a friend to help you.

Any hairs burned by the laser but not epilated can easily be removed using tweezers. Or let the hairs gradually fall out on their own. This may take up to 10 to 14 days. Red marks may appear after use, but these will gradually disappear within a few days following treatment.

The thicker the hair, the more effective the DLX-1088 at home Personal Laser Hair Remover

## 11. Q & A

### What is DLX-1088 laser hair removal?

Laser hair removal is the use of laser energy to produce permanent hair reduction. DLX-1088 emits light energy that is absorbed by the melanin pigment, resulting in rapid heating of the hair follicle. The heat destroys the papillae (hair follicles) that are responsible for hair growth.

### What areas can be treated with DLX-1088?

Any area, except adjacent to the eye, where there is excess hair can be treated with DLX-1088. The most common areas are the face, upper lip, neck, chest, periareolar, underarms, back, abdomen, bikini line, and legs.

### When will I see results from DLX-1088?

Most users see noticeable reduction after 4 treatments. It is important to understand the hair growth cycle. Hair grows in three different stages and only hairs in an active growth stage will be affected by a DLX-1088 treatment. This is one of the main reasons that multiple treatments are required to achieve the desired result.

### Are there any side effects?

In most cases there are no side effects associated with using DLX-1088. However, pigment change of the skin treated can rarely occur, including either dark or light areas. Pigment changes may persist for months, but are almost always temporary. They are very rare in patients with fair skin, who are un-tanned.

### Will laser hair removal work for me?

Due to the fact that a laser targets melanin, the more melanin an individual has in his/her hair, the more effective a laser will be. Light skin and dark hair are the best combination for laser hair removal, however DLX-1088 uses a diode laser which has a longer wavelength than other lasers used for hair removal, which may improve treatment results even on the darkest skin types.

These side effects include burns, blisters, and even hyper-pigmentation. But, the chance is less than 1% and it goes away soon. Extra care should be taken to avoid sun exposure after treatments with DLX-1088 for 2~3 days.

It is ok to take a shower after treatment. But avoid hot water in shower for 1 to 2 days. Don't treat the same area of skin more than once per hair removal session. Try to avoid overlapping treatments.

## 12. Troubleshooting

### 1. The device doesn't work when I put the skin contact bar to skins Check

That you plugged the power cord into the DLX-1088 and outlet. Unlock the laser with the safety key. Make sure Skin Contact Bar keeping contact with skins.

### 2. You hear the beep sound.

This beep sound that makes the go signal that the laser is now activated.

### 3. The device is getting hot.

The device can be getting hot by absorbing the heat if you use it for a long time.

### 13. Taking care of your DLX-1088 system

No special care is required. But it's a good idea to keep the DLX-1088 system components clean.

#### Cleaning

Unplug the adaptor from the DLX-1088 hand-piece and the power socket. And, Keep this device out of the reach of children. If needed, you should be able to remove dirt or debris easily with Q-tips. Be careful not to damage the lens of the laser.

#### Recycling/Disposal :

Please recycle or dispose of the devices properly.

Please follow the regulations where you live for the recycling or disposal of electrical devices.

Don't store the device in temperatures about 40 degree Celsius.(104 degree F)

### 14. Performance environment

Temperature :  $-10 < TC < 30\text{ }^{\circ}\text{C}$ ,

Relative Humidity :  $30\% < HR < 85\%$



SPECIFICATION			
Product	Laser Hair Remover (Class I Laser Product)	Model	DLX-1088
Wavelength	808nm	Size/Weight	188 x 60 x 68mm /142g
Material	ABS	Switching Power Adapter	220-240Vac 50Hz 0.3Amax DC6V 1A 4 x AA

